


Reficio 100[®] DBM Putty

DCI Donor Services Tissue Bank's Reficio 100[®] DBM Putty is comprised of 100% demineralized bone, maximizing allograft content and eliminating the need for extraneous carriers, fillers, or binding agents. The graft offers best in class handling and optimal regenerative capacity through a proprietary combination of demineralized cortical fibers and demineralized cortical particulate.

- Validated to produce a positive osteoinductive response¹
- Maintains moldable, cohesive handling characteristics to fill irregular defects and provide lavage resistance
- Engineered for enhanced regenerative capacity compared to particulate DBM²
- Sterility Assurance Level of 10⁻⁶ through low dose, low temperature gamma irradiation³
- Ready to use directly out of packaging with no need for graft preparation or hydration
- Proven regenerative capacity by exhibiting all five elements of new bone formation in a validated animal model⁴



Histological evidence of bone formation

DBM Putty	Size	Ambient
	1 cc	25800218
	2.5 cc	25812318
	5 cc	25800418
	10 cc	25800518



Reference:

1. Data on file, DCI Donor Services Tissue Bank
2. Martin GJ, et. Al. Spine, 1999, 24, 637-645
3. Data on file, DCI Donor Services Tissue Bank
4. Data on file, DCI Donor Services Tissue Bank

DCI Donor Services Tissue Bank is registered with the FDA, certified by the AATB and licensed with the states of California, Delaware, Florida, Illinois, Maryland, New York, and Oregon.

MK 109 Rev 001 2/2018

DCI Donor Services Tissue Bank
an extraordinary commitment to science, health and hope

